

Non elastic adhesive tape

- **DESCRIPTION :**

Non elastic strapping tape constructed from viscose fabric combined with a zinc oxide adhesive, Strappal® offers extra strong support and efficient stabilization of joint structures. This adhesive (free of natural rubber) helps to minimize the risk of allergic reaction.
- **CARACTERISTICS:**
 - Hand -tearable
 - Free of natural rubber.
 - Inextensible backing material provides rigid support
 - Serrated edges to prevent fraying
 - Sterilization not advisable
 - Consistent unwinding tension
- **INTENDED USES:**
 - Sport injuries (strains and sprains)
 - Orthopaedics (congenital dislocation of hip)
 - Skin adhesive traction.
 - Rigid immobilization of strained or sprained joints
 - Prophylaxis of sport injuries
- **DIRECTIONS FOR USE:**
 - To be applied to clean and dry skin
 - Recommended for patients sensitized to traditional adhesive tape containing natural rubber.
- **CAUTION FOR USE:**
 - Remove tape when the patient is submitted to radiological examination or radiation treatment (risk of opacity in X-ray negatives or diffraction of ionizing radiations in treatment).
- **CONTRA-INDICATIONS:**

Allergy to patient sensitized to synthetic products.

-
- **CE CLASS:**

Class I non sterile

First year of commercialization : 1983

- **COMPOSITION:**

Backing material Bleached viscose fabric
Average weight : 107 ± 7g/m²

Adhesives mass Composed of synthetic resin, lanolin, synthetic rubber, zinc oxide

Non elastic adhesive tape

Average weight: 105 ± 15 g/m²

- CONTROLS :

- of RAW MATERIALS:

Fabrics :

Weight-g/m², thread count, breaking load

Synthetic resin:

Solubility test , softening point, infra red spectroscopy.

Synthetic rubber

Solubility test, infra red spectroscopy

Lanoline: European Pharmacopoeia

Zinc oxide: density

- of FINISHED PRODUCT:

Adhesive mass weight

Peel adhesiveness control on steel plate

- STORAGE CONDITIONS: Keep away from light and heat (store below 25°C)
- SHELF LIFE: 3 years
- ADVISED METHOD OF DESTRUCTION :
Landfill according to local regulations.
- PACKAGING : Each strapping tape is rolled on cardboard core

Product code	Nart Code	Sizes	Packaging	Number per cardboard case
420100C	71496-00016-00	10m*4cm	1 / box	90
4202600	71488-00000-00	10m*2.5cm	8 by tube	160
4202650	71491-00000-00	10m*4cm	5 by tube	100
4202700	71494-00000-00	10m*5cm	4 by tube	80
4202755	71494-00007-00	10m*4cm	On tray	12
4202795	71496-00015-00	10m*1.25cm	On tray	36
4202800	71487-00000-00	10m*2.5cm	On tray	36
4202850	71490-00000-00	10m*4cm	On tray	24
4202900	71493-00000-00	10m*5cm	On tray	18
4202905	71494-00001-00	10m*7.5cm	On tray	6
48420050C	71496-00017-00	10m*2.5cm	1 / box	30
48420150	71489-00000-00	10m*2.5cm	1 / box Packed by 5	30
48420150C	71496-00018-00	10m*2.5cm	1 / box	30

Non elastic adhesive tape

Product code	Nart Code	Sizes	Presentation	Number per cardboard case
48420155	71492-00000-00	10m*4cm	1 / box Packed by 5	90
48420155C	71496-00019-00	10m*4cm	1 / box	90
72069-00011-00	72069-00011-00	5m*7.5cm	1 / box Packed by 12	96
72069-00012-00	72069-00012-00	5m*5cm	1 / box Packed by 12	144
72069-00013-00	72069-00013-00	5m*2.5cm	1 / box Packed by 12	216
71495-00001-00	71495-00001-00	10m*5xcm	1 / box	75

- LEGAL MANUFACTURER: BSN MEDICAL SAS
Rue du Millénaire –BP22
72320 Vibraye
FRANCE

Certificates : ISO 9001 - ISO 13485